

**RAPID TEST  
NOW AVAILABLE**

**ScheBo® • Biotech AG**  
a new era in life science

# ScheBo® • Pancreas Elastase1 Quick™

Stool Test for Exocrine Pancreatic Function Testing



## Faecal Pancreatic Elastase1

**The Gold Standard for Non-Invasive Exocrine Pancreatic Function Testing**

**... now also available as a rapid test**

The ScheBo® • Pancreas Elastase1 Quick™ stool test is a reliable and easy to perform rapid test for the detection of Pancreatic Elastase1 in stool samples.

Sample material: small stool sample

## Advantages of the

### ScheBo® • Pancreas Elastase1 Quick™ Test

- Test result available within minutes
- High sensitivity and specificity
- Small stool sample is sufficient
- Easy to perform
- No additional equipment necessary
- Pancreatic Elastase1 is stable during intestinal transit
- Absolutely pancreas-specific
- Uses monoclonal antibodies
- Enzyme therapy does not affect the test result
- Not affected by foodstuffs
- No special diet required

## Indications for use of the


### ScheBo® • Pancreas Elastase1 Quick™ Test

- Diagnosis/exclusion of pancreatic involvement in association with gastrointestinal symptoms: e.g., Maldigestion (Steatorrhoea/Diarrhoea), Abdominal Pain or 'Failure to Thrive'
- Diagnosis/exclusion of exocrine pancreatic insufficiency in: e.g., Chronic Pancreatitis, Cystic Fibrosis, Diabetes, Coeliac Disease, Irritable Bowel Syndrome (IBS), Inflammatory Bowel Diseases (IBD), Cholelithiasis (Gallstones), Osteoporosis or Papillary Stenosis

**Faecal Pancreatic Elastase1**  
**The Gold Standard for Non-Invasive**  
**Exocrine Pancreatic Function Testing**

**High Sensitivity and Specificity**

Diagnosis of Chronic Pancreatitis with the  
ScheBo® • Pancreatic Elastase1 Stool Test ELISA



Adapted from Löser, et al. Gut 1996;39:580-586

ScheBo® • Pancreas Elastase1 Quick™ has 92% sensitivity and 94% specificity when compared with the ScheBo® • Pancreatic Elastase1 Stool Test ELISA